

#### TESTIMONIAL

"It's great to see 21st Century technology that can replace 1940's technology. The McGRATH® MAC is really easy to become fluent with because it's just like a regular laryngoscope but has extra benefits over a direct laryngoscope. The enhanced on-screen view is a tremendous advantage while less force is exerted on the laryngeal structures. This device also creates the possibility of improving the grade of laryngoscopic view whilst using a conventional laryngoscopy technique."

STEPHEN LUNEY, M.D., CONSULTANT NEUROANESTHESIST, ROYAL VICTORIA HOSPITAL, BELFAST, UK

#### TESTIMONIAL

"For everyday practice it is a great device and has the potential to make us more successful every time."

IRENE OSBORN, M.D. DIRECTOR, DIVISION OF NEUROANESTHESIA, MOUNT SINAI MEDICAL CENTER, USA

#### TESTIMONIAL

"The McGRATH® MAC brings the most advanced airway management technology to routine intubation."

MICHAEL S. HIGGINS, PROFESSOR, DEPARTMENT OF ANESTHESIOLOGY, VANDERBILT UNIVERSITY MEDICAL CENTER

#### Technical Specifications

McGRATH®	MAC Handle
Size:	180mm x 68mm x 110mm
Weight:	200g
Power:	Proprietary 3.6V Lithium Battery Pack ( c. 250 minutes )
Light source:	High intensity LED
Display	2.5" LCD colour display
Camera:	CMOS
Materials:	Durable medical grade thermoplastics with reinforced structural alloy core. The device and packaging are latex free

#### McGRATH® MAC Disposable Blades

Material:	Medical grade optical polymer
Packaging:	Packaged sterile for single use

#### Regulatory

This product complies with EN 60601-1 and EN 60601-1-2 safety standards. The CE mark indicates that it meets the requirements of European Council Directives 93/42/EEC and 2007/47/EC concerning medical devices.

The device is regulated in the USA under FDA Regulation Number 868.5540 and device listed under the name McGRATH MAC.

"McGRATH" and "Aircraft" are registered trademarks of Aircraft Medical Limited. "CameraStick" is a trademark of Aircraft Medical Limited.

#### McGRATH® Video Laryngoscope range



McGRATH® Series 5 for specialists & difficult airways  
McGRATH® MAC for everyday use

#### Ordering

Contact your local seller for details (or visit [www.aircraftmedical.com](http://www.aircraftmedical.com))

AIRCRAFT  
MEDICAL



#### Aircraft Medical Limited

10 Saint Andrew Square, Edinburgh EH2 2AF United Kingdom  
tel: +44 (0) 131 718 6042 fax: +44 (0) 131 718 6100  
[www.aircraftmedical.com](http://www.aircraftmedical.com)  
COPYRIGHT 2011 © Aircraft Medical Limited Doc No. LE 004 v1.2

# McGRATH® MAC

The new video laryngoscope  
designed & priced for everyday use



# McGRATH® MAC

It's versatile. The McGRATH® MAC can be used as a direct or indirect view laryngoscope. It utilizes all new single-use McGRATH® MAC blades made from robust optical polymer. It's supported by an internal reinforced alloy chassis and CameraStick™ which delivers 'steel-like' rigidity. Its super-slim blade profile brings improved patient access and limits dental contact. Its unique, minute-by-minute countdown battery technology can last for weeks. It's compact, it's immersible and it's ready for the unexpected.

Video is here to stay, be part of it.



The highly refined McGRATH® MAC combines the strengths of the global standard Macintosh with McGRATH® technology.

**New McGRATH® MAC blade**  
Made from a robust optical polymer and supported by an internal reinforced CameraStick™, the McGRATH® MAC slimline single-use blade delivers a 'steel-like' rigidity.

