

# *L770 Bag Valve Masks*




## *The Evolution of Resuscitation Continues.*

The LSP L770 series of disposable bag valve masks offers accurate dependable operation coupled with the ability to manually feel a patient's respiratory response from a safe distance.

By separating the bag from the patient valve, the forces created by squeezing the bag are not transferred to the mask or ET tube. This helps to keep a tight seal when using a mask, and prevents dislodging the ET tube when bagging intubated patients. The separation also allows the person bagging to walk alongside the patient during transport, ensuring more consistent delivery of tidal volumes to the patient and allowing the transporters to move quickly.

In an emergency vehicle, the caregiver can remain safely seated and strapped in during transit. Using a remote bag creates more space around the patient's head, yet you can still see the chest rise and fall. The added tubing does not affect delivered tidal volume or oxygen concentration.

# L770 Bag Valve Masks



## ORDERING INFORMATION

- L770-040** Adult Bag Valve Mask with reservoir, mask, flex tubing, 7' oxygen tubing (no pressure relief)
- L770-100** Child Bag Valve Mask with reservoir, mask, flex tubing, 7' oxygen tubing, pressure relief
- L770-101** Child Bag Valve Mask with reservoir, mask, flex tubing, 7' oxygen tubing (no pressure relief)
- L770-200** Infant Bag Valve Mask with reservoir, mask, flex tubing, 7' oxygen tubing, pressure relief
- L770-202** Infant/Neonate Bag Valve Mask with reservoir, 2 masks, flex tubing, 7' oxygen tubing, pressure relief

## ACCESSORIES

- 890113** Disposable Mask Restraint System (12/case)
- 890060** Disposable Manometer (12/case)
- LPEEP** Disposable PEEP Valve with adapter (12/case)



Allied Healthcare Products, Inc. • 1720 Sublette Avenue • St. Louis, MO 63110 USA  
314/771-2400 • 800/444-3954 • FAX: 314/771-0650 • [www.alliedhpi.com](http://www.alliedhpi.com)

U.S.  
Toll Free: 800/444-3940  
FAX: 800/477-7701

Canada & Mexico  
Toll Free: 800/446-0552  
FAX: 800/246-6201

International  
Telephone: 314/268-1683  
FAX: 314/771-5183

**ISO 13485: 2003**

©Allied Healthcare Products, Inc.®

All specifications are nominal and  
subject to change without notice.

Form 90-00-0006  
Rev. July 2010