






8 SPRING STREET, STAFFORD SPRINGS, CT 06076 • 800 344 6227 • WWW.HOBBSMEDICAL.COM



HOBBS MEDICAL, INC.



INTRODUCTION

WELCOME

Established in 1984, Hobbs Medical, Inc.

has built a solid reputation for

innovation and quality. A privately

held company with expertise in the

development and manufacture of

diagnostic and therapeutic accessories for

flexible GI and Pulmonary Endoscopy, Hobbs

works with a network of industry professionals,

consultants, sales representatives and international distribution agencies

to bring precision endoscopic accessories to the world market. We take pride in

our comprehensive line of precision accessories, all manufactured in the USA, and

continue to introduce unique devices to the ever-changing market.

The entire staff at Hobbs Medical is committed to exceeding each and every customer's expectations;

providing cost-effective, quality and innovative products and investing in ongoing product development.

Our Sales personnel are always available to answer inquiries. We look forward to responding to your needs.

Contact Sales at (800) 344-6227 or visit www.hobbsmedical.com



HOBBS MEDICAL, INC. 8 SPRING STREET, STAFFORD SPRINGS, CT 06076 • 800 344 6227 • WWW.HOBBSMEDICAL.COM

ALL HOBBS MEDICAL PRODUCTS ARE PROUDLY MANUFACTURED IN THE USA.
HOBBS MEDICAL IS ISO-13485 REGISTERED; CE-CERTIFIED AND COMPLIES WITH
COUNCIL DIRECTIVE 93/94EEC AND FDA QUALITY SYSTEMS REGULATIONS.

Customer Service
800-344-6227 toll-free
860-684-5875
860-684-7574 fax
customerservice@hobbsmedical.com

Sales
800-344-6227 toll-free
860-684-5875
salesdept@hobbsmedical.com

QUALITY POLICY

The Quality Policy of Hobbs Medical, Inc. is to meet the needs and expectations of our customers. This will be accomplished through market feedback and continuous evaluation of our Quality Management System to ensure its continued effectiveness and suitability in complying with customer needs and all regulatory requirements at all levels of the organization.

SPECIFICATION CHANGES

Hobbs Medical reserves the right to make improvements and changes in the products described herein. Should you have a suggestion for a new device or an improvement of an existing device, let our Engineering staff translate your idea into reality. Our website (www.hobbsmedical.com) provides an avenue for product suggestions.

PRODUCT WARRANTY

Descriptions or specifications appearing in Hobbs Medical literature are solely meant to generally describe the product at the time of manufacture and do not constitute any express warranties. Hobbs Medical warrants all of its products to be free from defects in materials and workmanship. Storage conditions and product use are factors beyond the control of Hobbs Medical and may directly affect the products and the results obtained from them. The user assumes responsibility for any consequences, which may arise from the final use of the products. No person has any authority to bind Hobbs Medical to any representation or warranty concerning these products.

All reusable products in this catalog are warranted to the original purchaser to be free from defects in material and workmanship for a period of six months from the date of purchase. All disposable products are warranted to the original purchaser to be free from defects in materials and workmanship for the first use only.

RETURN POLICY

A Return Authorization Number must be obtained for all returned products. Contact our Customer Service Department (800-344-6227). Returns must include a reason for return, lot number, catalog number and invoice number. Contaminated products sent for evaluation will not be accepted unless the product has been cleaned and high level disinfected prior to returning to Hobbs Medical. Returned, unused products must be in the original unopened Hobbs Medical pouch and will be subject to a 20% restocking fee.

SHIPPING & TERMS

Most items are shipped within one day. Products are shipped via ground transportation unless otherwise specified. MasterCard and VISA accepted. F.O.B. Stafford Springs, Connecticut. Due to fluctuating operating costs, we reserve the right to change prices without prior notice.

ORDERING INFORMATION

Our Sales and Customer Service personnel are always available to answer inquiries and to meet our customers' specific needs. Order and Price Requests can be placed by telephone, fax or email.

CONTACT INFORMATION

Customer Service
800-344-6227 toll-free
860-684-5875
860-684-7574 fax
customerservice@hobbsmedical.com

Sales
800-344-6227 toll-free
860-684-5875
salesdept@hobbsmedical.com

Orders & Returns

Product Information
Pricing
Specification Changes

BIOPSY FORCEPS

BIOPSY FORCEPS

GASTROINTESTINAL PRODUCTS

Cold Biopsy Forceps: Reusable

- A variety of jaw configurations and lengths to meet your procedural requirements
- Needle forceps provide several biopsies with one-pass for efficient specimen retrieval
- Precision milled biopsy cups ensure clean-cut, non-deformed specimens
- Smoothly polished tapered distal end reduces channel friction during procedural maneuvers
- Validated for reprocessing – see instructions for use
- New Endoscope compatible lengths
- No minimum order requirements
- Packaged non-sterile



COLD BIOPSY FORCEPS: REUSABLE

Catalog Number	Description	Cup	Sheath	Length
12404	Fenestrated	Oval	1.8mm	160cm
12448	Needle	Oval	1.8mm	160cm
12440	Fenestrated	Oval	2.3mm	100cm
12442	Fenestrated	Oval	2.3mm	180cm
12450	Needle	Oval	2.3mm	180cm
12825	Serrated	Alligator	2.3mm	180cm
12444	Fenestrated	Oval	2.3mm	230cm
12452	Needle	Oval	2.3mm	230cm
12827	Serrated	Alligator	2.3mm	230cm
12446	Fenestrated	Oval	2.3mm	320cm
12454	Needle	Oval	2.3mm	320cm

PULMONARY PRODUCTS

COLD BIOPSY FORCEPS: REUSABLE


Catalog Number	Description	Cup	Sheath	Length
12402	Fenestrated	Oval	1.8mm	120cm
12821	Serrated	Alligator	1.8mm	120cm



GASTROINTESTINAL PRODUCTS

Hot Biopsy Forceps: Reusable

- Precision milled biopsy cups ensure clean-cut, non-deformed cauterized specimens
- Handle features a rotatable active cord fitting
- Recommended for a 2.8mm or larger working channel
- Validated for reprocessing – see instructions for use
- Requires Active Cord*
- No minimum order requirements
- Packaged non-sterile



*See Specialty Items, page 34 for Active Cord

HOT BIOPSY FORCEPS: REUSABLE

Catalog Number	Description	Cup	Sheath	Length
13486NH	No Handle	Oval	2.3mm	180cm
13486WH	With Handle	Oval	2.3mm	180cm
13487NH	No Handle	Oval	2.3mm	230cm
13487WH	With Handle	Oval	2.3mm	230cm

HOT BIOPSY FORCEPS REPLACEMENT HANDLE: REUSABLE

Catalog Number	Description
4506B	Hot Biopsy Handle

