




HOBBS MEDICAL, INC.

8 SPRING STREET, STAFFORD SPRINGS, CT 06076 • 800 344 6227 • WWW.HOBBSMEDICAL.COM



HOBBS MEDICAL, INC.



## INTRODUCTION

### WELCOME

Established in 1984, Hobbs Medical, Inc.

has built a solid reputation for

innovation and quality. A privately

held company with expertise in the

development and manufacture of

diagnostic and therapeutic accessories for

flexible GI and Pulmonary Endoscopy, Hobbs

works with a network of industry professionals,

consultants, sales representatives and international distribution agencies

to bring precision endoscopic accessories to the world market. We take pride in

our comprehensive line of precision accessories, all manufactured in the USA, and

continue to introduce unique devices to the ever-changing market.

The entire staff at Hobbs Medical is committed to exceeding each and every customer's expectations;

providing cost-effective, quality and innovative products and investing in ongoing product development.

Our Sales personnel are always available to answer inquiries. We look forward to responding to your needs.

Contact Sales at (800) 344-6227 or visit [www.hobbsmedical.com](http://www.hobbsmedical.com)



**HOBBS MEDICAL, INC.** 8 SPRING STREET, STAFFORD SPRINGS, CT 06076 • 800 344 6227 • [WWW.HOBBSMEDICAL.COM](http://WWW.HOBBSMEDICAL.COM)



ALL HOBBS MEDICAL PRODUCTS ARE PROUDLY MANUFACTURED IN THE USA.  
HOBBS MEDICAL IS ISO-13485 REGISTERED; CE-CERTIFIED AND COMPLIES WITH  
COUNCIL DIRECTIVE 93/94EEC AND FDA QUALITY SYSTEMS REGULATIONS.

**Customer Service**  
800-344-6227 toll-free  
860-684-5875  
860-684-7574 fax  
[customerservice@hobbsmedical.com](mailto:customerservice@hobbsmedical.com)

**Sales**  
800-344-6227 toll-free  
860-684-5875  
[salesdept@hobbsmedical.com](mailto:salesdept@hobbsmedical.com)

#### QUALITY POLICY

The Quality Policy of Hobbs Medical, Inc. is to meet the needs and expectations of our customers. This will be accomplished through market feedback and continuous evaluation of our Quality Management System to ensure its continued effectiveness and suitability in complying with customer needs and all regulatory requirements at all levels of the organization.

#### SPECIFICATION CHANGES

Hobbs Medical reserves the right to make improvements and changes in the products described herein. Should you have a suggestion for a new device or an improvement of an existing device, let our Engineering staff translate your idea into reality. Our website ([www.hobbsmedical.com](http://www.hobbsmedical.com)) provides an avenue for product suggestions.

#### PRODUCT WARRANTY

Descriptions or specifications appearing in Hobbs Medical literature are solely meant to generally describe the product at the time of manufacture and do not constitute any express warranties. Hobbs Medical warrants all of its products to be free from defects in materials and workmanship. Storage conditions and product use are factors beyond the control of Hobbs Medical and may directly affect the products and the results obtained from them. The user assumes responsibility for any consequences, which may arise from the final use of the products. No person has any authority to bind Hobbs Medical to any representation or warranty concerning these products.

All reusable products in this catalog are warranted to the original purchaser to be free from defects in material and workmanship for a period of six months from the date of purchase. All disposable products are warranted to the original purchaser to be free from defects in materials and workmanship for the first use only.

#### RETURN POLICY

A Return Authorization Number must be obtained for all returned products. Contact our Customer Service Department (800-344-6227). Returns must include a reason for return, lot number, catalog number and invoice number. Contaminated products sent for evaluation will not be accepted unless the product has been cleaned and high level disinfected prior to returning to Hobbs Medical. Returned, unused products must be in the original unopened Hobbs Medical pouch and will be subject to a 20% restocking fee.

#### SHIPPING & TERMS

Most items are shipped within one day. Products are shipped via ground transportation unless otherwise specified. MasterCard and VISA accepted. F.O.B. Stafford Springs, Connecticut. Due to fluctuating operating costs, we reserve the right to change prices without prior notice.

#### ORDERING INFORMATION

Our Sales and Customer Service personnel are always available to answer inquiries and to meet our customers' specific needs. Order and Price Requests can be placed by telephone, fax or email.

#### CONTACT INFORMATION

**Customer Service**  
800-344-6227 toll-free  
860-684-5875  
860-684-7574 fax  
[customerservice@hobbsmedical.com](mailto:customerservice@hobbsmedical.com)

*Orders & Returns*

**Sales**  
800-344-6227 toll-free  
860-684-5875  
[salesdept@hobbsmedical.com](mailto:salesdept@hobbsmedical.com)

*Product Information  
Pricing  
Specification Changes*

# POLYPECTOMY SNARES

POLYPECTOMY SNARES

## GASTROINTESTINAL PRODUCTS

### Polypectomy Snares: Reusable

- Handle features a rotatable Active Cord fitting\*
- Fully rotatable handle for ease of manipulation and placement
- Recommended for a 2.8mm or larger working channel
- Validated for reprocessing – see instructions for use
- Requires Active Cord\*
- No minimum order requirements



\*See Specialty Items, page 34 for Active Cord

#### POLYPECTOMY SNARES: REUSABLE (Snare and sheath assembly only – handle sold separately)

Catalog Number	Description	Loop Diameter	Sheath	Length	Connector
4510	Oval Loop	4.0cm x 3.0cm	2.3mm	220cm	Female
4511	Oval Loop	5.0cm x 3.5cm	2.3mm	220cm	Female
4514	Hexagonal Loop	4.0cm x 3.5cm	2.3mm	220cm	Female
4518	Crescent Loop	3.5cm x 3.5cm	2.3mm	220cm	Female

#### ROTATABLE POLYPECTOMY REPLACEMENT HANDLE: REUSABLE

Catalog Number	Description
4506	Snare Handle

#### REUSABLE POLYPECTOMY SNARE KITS

Catalog Number	Description	Loop Diameter	Sheath	Length	Connector
4501	Oval Loop	4.0cm x 3.0cm	2.3mm	220cm	Female
4503	Hexagonal Loop	4.0cm x 3.5cm	2.3mm	220cm	Female
4505	Crescent Loop	3.5cm x 3.5cm	2.3mm	220cm	Female
4507	Oval Loop	5.0cm x 3.5cm	2.3mm	220cm	Female

#### Includes:


- Two Snares
- One Handle – #4506
- One Active Cord – Identify brand of cautery unit in order to include the appropriate cord



## GASTROINTESTINAL PRODUCTS

### Polypectomy Snares: Disposable

- Braided stainless steel wire prevents loop distortion
- Both male and female cautery connections offered to meet individual cautery needs
- Completely disposable unit prevents cross-contamination concerns
- Requires Active Cord\*
- Packaged sterile: 5 per box




\*See Specialty Items, page 34 for Active Cord

#### POLYPECTOMY SNARES: DISPOSABLE

Catalog Number	Description	Loop Diameter	Sheath	Length	Connector
4561	Hexagonal Loop	4.5cm x 3.0cm	2.3mm	220cm	Female
4561-OLY	Hexagonal Loop	4.5cm x 3.0cm	2.3mm	220cm	Olympus
4562	Crescent Loop	5.0cm x 2.5cm	2.3mm	220cm	Female
4562-OLY	Crescent Loop	5.0cm x 2.5cm	2.3mm	220cm	Olympus
4563	Standard Oval Loop	5.0cm x 2.5cm	2.3mm	220cm	Female
4563-OLY	Standard Oval Loop	5.0cm x 2.5cm	2.3mm	220cm	Olympus
4564	Maxi Oval Loop	5.5cm x 3.5cm	2.3mm	220cm	Female
4564-OLY	Maxi Oval Loop	5.5cm x 3.5cm	2.3mm	220cm	Olympus
4565	Mini Oval Loop	3.5cm x 2.0cm	2.3mm	220cm	Female
4565-OLY	Mini Oval Loop	3.5cm x 2.0cm	2.3mm	220cm	Olympus
4567	Micro Oval Loop	2.5cm x 1.5cm	2.3mm	220cm	Female
4567-OLY	Micro Oval Loop	2.5cm x 1.5cm	2.3mm	220cm	Olympus



## GASTROINTESTINAL PRODUCTS



### Polypectomy Snares: Short Throw Disposable

- Handle designed for comfortable extension during polypectomy procedures
- Completely disposable unit prevents cross-contamination concerns
- Braided stainless steel wire prevents loop distortion
- Both male and female cautery connections offered to meet individual cautery needs
- Standard Oval available for a 2.0mm working channel
- Requires Active Cord\*
- Packaged sterile: 5 per box

\*See Specialty Items, page 34 for Active Cord

#### POLYPECTOMY SNARES: SHORT THROW DISPOSABLE

Catalog Number	Description	Loop Diameter	Sheath	Length	Connector
4571	Hexagonal Loop	4.5cm x 3.0cm	2.3mm	220cm	Female
4571-OLY	Hexagonal Loop	4.5cm x 3.0cm	2.3mm	220cm	Olympus
4572	Crescent Loop	5.0cm x 2.5cm	2.3mm	220cm	Female
4572-OLY	Crescent Loop	5.0cm x 2.5cm	2.3mm	220cm	Olympus
4573	Standard Oval Loop	5.0cm x 2.5cm	2.3mm	220cm	Female
4573-OLY	Standard Oval Loop	5.0cm x 2.5cm	2.3mm	220cm	Olympus
4575	Mini Oval Loop	3.5cm x 2.0cm	2.3mm	220cm	Female
4575-OLY	Mini Oval Loop	3.5cm x 2.0cm	2.3mm	220cm	Olympus
4577	Micro Oval Loop	2.5cm x 1.5cm	2.3mm	220cm	Female
4577-OLY	Micro Oval Loop	2.5cm x 1.5cm	2.3mm	220cm	Olympus
7202	Oval Loop	3.5cm x 2.0cm	1.8mm	220cm	Female

#### CONTACT INFORMATION

**Customer Service**  
 800-344-6227 toll-free  
 860-684-5875  
 860-684-7574 fax  
 customerservice@hobbsmedical.com

*Orders & Returns*

**Sales**  
 800-344-6227 toll-free  
 860-684-5875  
 salesdept@hobbsmedical.com

*Product Information  
 Pricing  
 Specification Changes*

