






8 SPRING STREET, STAFFORD SPRINGS, CT 06076 • 800 344 6227 • WWW.HOBBSMEDICAL.COM



HOBBS MEDICAL, INC.



INTRODUCTION

WELCOME

Established in 1984, Hobbs Medical, Inc.

has built a solid reputation for

innovation and quality. A privately

held company with expertise in the

development and manufacture of

diagnostic and therapeutic accessories for

flexible GI and Pulmonary Endoscopy, Hobbs

works with a network of industry professionals,

consultants, sales representatives and international distribution agencies

to bring precision endoscopic accessories to the world market. We take pride in

our comprehensive line of precision accessories, all manufactured in the USA, and

continue to introduce unique devices to the ever-changing market.

The entire staff at Hobbs Medical is committed to exceeding each and every customer's expectations;

providing cost-effective, quality and innovative products and investing in ongoing product development.

Our Sales personnel are always available to answer inquiries. We look forward to responding to your needs.

Contact Sales at (800) 344-6227 or visit www.hobbsmedical.com



HOBBS MEDICAL, INC.

8 SPRING STREET, STAFFORD SPRINGS, CT 06076 • 800 344 6227 • WWW.HOBBSMEDICAL.COM

ALL HOBBS MEDICAL PRODUCTS ARE PROUDLY MANUFACTURED IN THE USA.
HOBBS MEDICAL IS ISO-13485 REGISTERED; CE-CERTIFIED AND COMPLIES WITH
COUNCIL DIRECTIVE 93/94EEC AND FDA QUALITY SYSTEMS REGULATIONS.

Customer Service
800-344-6227 toll-free
860-684-5875
860-684-7574 fax
customerservice@hobbsmedical.com

Sales
800-344-6227 toll-free
860-684-5875
salesdept@hobbsmedical.com

QUALITY POLICY

The Quality Policy of Hobbs Medical, Inc. is to meet the needs and expectations of our customers. This will be accomplished through market feedback and continuous evaluation of our Quality Management System to ensure its continued effectiveness and suitability in complying with customer needs and all regulatory requirements at all levels of the organization.

SPECIFICATION CHANGES

Hobbs Medical reserves the right to make improvements and changes in the products described herein. Should you have a suggestion for a new device or an improvement of an existing device, let our Engineering staff translate your idea into reality. Our website (www.hobbsmedical.com) provides an avenue for product suggestions.

PRODUCT WARRANTY

Descriptions or specifications appearing in Hobbs Medical literature are solely meant to generally describe the product at the time of manufacture and do not constitute any express warranties. Hobbs Medical warrants all of its products to be free from defects in materials and workmanship. Storage conditions and product use are factors beyond the control of Hobbs Medical and may directly affect the products and the results obtained from them. The user assumes responsibility for any consequences, which may arise from the final use of the products. No person has any authority to bind Hobbs Medical to any representation or warranty concerning these products.

All reusable products in this catalog are warranted to the original purchaser to be free from defects in material and workmanship for a period of six months from the date of purchase. All disposable products are warranted to the original purchaser to be free from defects in materials and workmanship for the first use only.

RETURN POLICY

A Return Authorization Number must be obtained for all returned products. Contact our Customer Service Department (800-344-6227). Returns must include a reason for return, lot number, catalog number and invoice number. Contaminated products sent for evaluation will not be accepted unless the product has been cleaned and high level disinfected prior to returning to Hobbs Medical. Returned, unused products must be in the original unopened Hobbs Medical pouch and will be subject to a 20% restocking fee.

SHIPPING & TERMS

Most items are shipped within one day. Products are shipped via ground transportation unless otherwise specified. MasterCard and VISA accepted. F.O.B. Stafford Springs, Connecticut. Due to fluctuating operating costs, we reserve the right to change prices without prior notice.

ORDERING INFORMATION

Our Sales and Customer Service personnel are always available to answer inquiries and to meet our customers' specific needs. Order and Price Requests can be placed by telephone, fax or email.


CONTACT INFORMATION

Customer Service
800-344-6227 toll-free
860-684-5875
860-684-7574 fax
customerservice@hobbsmedical.com

Sales
800-344-6227 toll-free
860-684-5875
salesdept@hobbsmedical.com

Orders & Returns

Product Information
Pricing
Specification Changes



GASTROINTESTINAL PRODUCTS

**Microbiology Brushes:
Protected Specimen Brush**

- Teflon sheathing for ease of insertion and manipulation during procedures
- Spring-loaded handle and permanent tear drop shaped tip ensures a protective seal
- Latch hook for locking handle in open brush position
- Ergonomic handle provides operator comfort
- Brush configuration ensures a high specimen yield
- Packaged sterile for single use: 5 per box



MICROBIOLOGY BRUSHES

Catalog Number	Brush Diameter	Sheath	Length	Handle Color
4352	2.0mm	1.7mm	220cm	Green Swirl

PULMONARY PRODUCTS

MICROBIOLOGY BRUSHES


Catalog Number	Brush Diameter	Sheath	Length	Handle Color
4310	1.0mm	1.7mm	100cm	Blue Swirl
4320	2.0mm	1.7mm	100cm	Blue Swirl



GASTROINTESTINAL AND PULMONARY PRODUCTS

**Modified Saccomanno™
Cell Transport System**

- Allows for quick washing of brushes, forceps and needles
- Immediate cell fixation
- Convenient self-contained vial with patient ID tags
- Packaged 10 per box



MODIFIED SACCOMANNO™ CELL TRANSPORT SYSTEM


Catalog Number	Description
4101	Cell Transport System



GASTROINTESTINAL PRODUCTS

Cytology Brushes

- Teflon sheathing for ease of insertion and manipulation during procedures
- High tempered stainless steel wire to prevent kinking
- Radiopaque ball tip prevents mucosal trauma
- Thumb ring for easy operation
- Endoscope compatible lengths
- Packaged sterile for single use: 10 per box



CYTOLOGY BRUSHES

Catalog Number	Brush Diameter	Sheath	Length	Handle Color
4102	3.0mm	1.7mm	160cm	Green
4104	3.0mm	1.7mm	200cm	Green
4106	3.0mm	1.7mm	230cm	Green
4108	3.0mm	2.3mm	160cm	Green
4110	3.0mm	2.3mm	200cm	Green
4112	3.0mm	2.3mm	260cm	Green

PULMONARY PRODUCTS

CYTOLOGY BRUSHES

Catalog Number	Brush Diameter	Sheath	Length	Handle Color
4202	1.0mm	1.7mm	120cm	Blue
4204	2.0mm	1.7mm	120cm	Blue
4206	3.0mm	1.7mm	120cm	Blue

NASAL CYTOLOGY BRUSH

Catalog Number	Brush Diameter	Sheath	Length	Handle Color
4290	2.0mm	N/A	10cm	Blue Swirl




GASTROINTESTINAL PRODUCTS

Cleaning Brushes: Disposable

GOOD-BI-BRUSH (2-BRUSH CLEANING SYSTEM)

- Unique tapered 2-brush system provides dual cleaning action with one pass
- Blunt, looped tip brush protects delicate channels
- Improved sheath flexibility for easy passage through working channel and umbilical exit port
- Vivid color cannula offers easy identification in cleaning solution
- Endoscope compatible lengths
- Packaged for single use: 25 per box



GOOD-BI-BRUSH (2-BRUSH CLEANING SYSTEM): DISPOSABLE

Catalog Number	Brush Taper	Brush Length	Sheath	Length	Working Channel
4660	3-5mm x 4-6mm	2cm x 1cm	1.7mm	260cm	2.0mm & Larger

Single Cleaning Brush System

- Unique tapered brush provides greater channel cleaning

SINGLE CLEANING BRUSH SYSTEM: DISPOSABLE

Catalog Number	Brush Taper	Brush Length	Sheath	Length	Working Channel
4676	3-5mm	2cm	1.7mm	260cm	2.0mm & Larger
4677	3-5mm	1cm	1.7mm	260cm	2.0mm & Larger

PULMONARY PRODUCTS

- Unique tapered 2-brush system provides dual cleaning action with one pass

GOOD-BI-BRUSH (2-BRUSH CLEANING SYSTEM): DISPOSABLE

Catalog Number	Brush Taper	Brush Length	Sheath	Length	Working Channel
4665	3-5mm x 4-6mm	2cm x 1cm	1.7mm	100cm	2.0mm & Larger

