






8 SPRING STREET, STAFFORD SPRINGS, CT 06076 • 800 344 6227 • WWW.HOBBSMEDICAL.COM



HOBBS MEDICAL, INC.



INTRODUCTION

WELCOME

Established in 1984, Hobbs Medical, Inc.

has built a solid reputation for

innovation and quality. A privately

held company with expertise in the

development and manufacture of

diagnostic and therapeutic accessories for

flexible GI and Pulmonary Endoscopy, Hobbs

works with a network of industry professionals,

consultants, sales representatives and international distribution agencies

to bring precision endoscopic accessories to the world market. We take pride in

our comprehensive line of precision accessories, all manufactured in the USA, and

continue to introduce unique devices to the ever-changing market.

The entire staff at Hobbs Medical is committed to exceeding each and every customer's expectations;

providing cost-effective, quality and innovative products and investing in ongoing product development.

Our Sales personnel are always available to answer inquiries. We look forward to responding to your needs.

Contact Sales at (800) 344-6227 or visit www.hobbsmedical.com



HOBBS MEDICAL, INC. 8 SPRING STREET, STAFFORD SPRINGS, CT 06076 • 800 344 6227 • WWW.HOBBSMEDICAL.COM

ALL HOBBS MEDICAL PRODUCTS ARE PROUDLY MANUFACTURED IN THE USA.
HOBBS MEDICAL IS ISO-13485 REGISTERED; CE-CERTIFIED AND COMPLIES WITH
COUNCIL DIRECTIVE 93/94EEC AND FDA QUALITY SYSTEMS REGULATIONS.

Customer Service
800-344-6227 toll-free
860-684-5875
860-684-7574 fax
customerservice@hobbsmedical.com

Sales
800-344-6227 toll-free
860-684-5875
salesdept@hobbsmedical.com

QUALITY POLICY

The Quality Policy of Hobbs Medical, Inc. is to meet the needs and expectations of our customers. This will be accomplished through market feedback and continuous evaluation of our Quality Management System to ensure its continued effectiveness and suitability in complying with customer needs and all regulatory requirements at all levels of the organization.

SPECIFICATION CHANGES

Hobbs Medical reserves the right to make improvements and changes in the products described herein. Should you have a suggestion for a new device or an improvement of an existing device, let our Engineering staff translate your idea into reality. Our website (www.hobbsmedical.com) provides an avenue for product suggestions.

PRODUCT WARRANTY

Descriptions or specifications appearing in Hobbs Medical literature are solely meant to generally describe the product at the time of manufacture and do not constitute any express warranties. Hobbs Medical warrants all of its products to be free from defects in materials and workmanship. Storage conditions and product use are factors beyond the control of Hobbs Medical and may directly affect the products and the results obtained from them. The user assumes responsibility for any consequences, which may arise from the final use of the products. No person has any authority to bind Hobbs Medical to any representation or warranty concerning these products.

All reusable products in this catalog are warranted to the original purchaser to be free from defects in material and workmanship for a period of six months from the date of purchase. All disposable products are warranted to the original purchaser to be free from defects in materials and workmanship for the first use only.

RETURN POLICY

A Return Authorization Number must be obtained for all returned products. Contact our Customer Service Department (800-344-6227). Returns must include a reason for return, lot number, catalog number and invoice number. Contaminated products sent for evaluation will not be accepted unless the product has been cleaned and high level disinfected prior to returning to Hobbs Medical. Returned, unused products must be in the original unopened Hobbs Medical pouch and will be subject to a 20% restocking fee.

SHIPPING & TERMS

Most items are shipped within one day. Products are shipped via ground transportation unless otherwise specified. MasterCard and VISA accepted. F.O.B. Stafford Springs, Connecticut. Due to fluctuating operating costs, we reserve the right to change prices without prior notice.

ORDERING INFORMATION

Our Sales and Customer Service personnel are always available to answer inquiries and to meet our customers' specific needs. Order and Price Requests can be placed by telephone, fax or email.

CONTACT INFORMATION

Customer Service
800-344-6227 toll-free
860-684-5875
860-684-7574 fax
customerservice@hobbsmedical.com

Sales
800-344-6227 toll-free
860-684-5875
salesdept@hobbsmedical.com

Orders & Returns

Product Information
Pricing
Specification Changes


RETRIEVAL DEVICES

RETRIEVAL DEVICES

GASTROINTESTINAL PRODUCTS

Helical Retrieval Baskets: Reusable

- 3-ring handle designed for comfortable extension during procedures
- Luer-lock side port for fluid injection
- Compatible with 2.8mm and larger channels
- Choice of open basket diameter or wire configuration
- Replacement handles available
- Validated for reprocessing – see instructions for use
- No minimum order requirements
- Packaged sterile



HELICAL RETRIEVAL BASKETS: REUSABLE (Packaged fully assembled with handle)

Catalog Number	Description	Open Diameter	Sheath	Length
4830	4-Wire	15mm	2.3mm	220cm
4832	4-Wire	20mm	2.3mm	220cm
4814	3-Wire	15mm	2.3mm	220cm
4816	3-Wire	20mm	2.3mm	220cm

HELICAL RETRIEVAL BASKET REPLACEMENT HANDLE: REUSABLE

Catalog Number	Description
4850	Basket Handle

PULMONARY PRODUCTS

- Compatible with 2.0mm and larger channels

HELICAL RETRIEVAL BASKETS: REUSABLE


Catalog Number	Description	Open Diameter	Sheath	Length
4806	3-Wire	15mm	1.8mm	150cm
4808	3-Wire	20mm	1.8mm	150cm
4822	4-Wire	15mm	1.8mm	150cm
4824	4-Wire	20mm	1.8mm	150cm



GASTROINTESTINAL PRODUCTS

Retrieval Forceps: Reusable

- 3-ring handle designed for comfortable extension during procedures
- Luer-lock side port for fluid injection
- Replacement handles available
- Validated for reprocessing – see instructions for use
- No minimum order requirements
- Compatible with 2.8mm and larger channels
- Packaged sterile



RETRIEVAL FORCEPS: REUSABLE (Packaged fully assembled with handle)

Catalog Number	Description	Open Diameter	Sheath	Length
5008	3-Prong	20mm	2.3mm	220cm
5018	4-Prong	20mm	2.3mm	220cm

RETRIEVAL FORCEPS REPLACEMENT HANDLE: REUSABLE

Catalog Number	Description
5050	Forceps Handle

PULMONARY PRODUCTS

- Compatible with 2.0mm and larger channels

RETRIEVAL FORCEPS: REUSABLE


Catalog Number	Description	Open Diameter	Sheath	Length
5004	3-Prong	20mm	1.8mm	150cm



GASTROINTESTINAL PRODUCTS

**Helical Retrieval Baskets:
Reusable**

- 3-ring handle designed for comfortable extension during procedures
- Luer-lock side port for fluid injection
- Compatible with 2.8mm and larger channels
- Choice of open basket diameter or wire configuration
- Replacement handles available
- Validated for reprocessing – see instructions for use
- No minimum order requirements
- Packaged sterile



HELICAL RETRIEVAL BASKETS: REUSABLE (Packaged fully assembled with handle)

Catalog Number	Description	Open Diameter	Sheath	Length
4830	4-Wire	15mm	2.3mm	220cm
4832	4-Wire	20mm	2.3mm	220cm
4814	3-Wire	15mm	2.3mm	220cm
4816	3-Wire	20mm	2.3mm	220cm

HELICAL RETRIEVAL BASKET REPLACEMENT HANDLE: REUSABLE

Catalog Number	Description
4850	Basket Handle

PULMONARY PRODUCTS

- Compatible with 2.0mm and larger channels

HELICAL RETRIEVAL BASKETS: REUSABLE

Catalog Number	Description	Open Diameter	Sheath	Length
4806	3-Wire	15mm	1.8mm	150cm
4808	3-Wire	20mm	1.8mm	150cm
4822	4-Wire	15mm	1.8mm	150cm
4824	4-Wire	20mm	1.8mm	150cm



CONTACT INFORMATION

Customer Service
800-344-6227 toll-free
860-684-5875
860-684-7574 fax
customerservice@hobbsmedical.com

Orders & Returns

Sales
800-344-6227 toll-free
860-684-5875
salesdept@hobbsmedical.com

*Product Information
Pricing
Specification Changes*